Instructions for Use

 Manufacturer: C-Pro Direct Ltd

 Catalogue No: DORSIRAMP

Black Robin Designs

C-Pro Direct Ltd

Unit 2, Black Robins Farm,

Grants Lane, Edenbridge, Kent, TN8 6QP, UK.

Tel: +44 (0) 1732 860 158 www.c-prodirect.com

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GB Instructions for Use

Intended Use

Dorsi Ramp is a physiotherapy device to help children to stretch their Achilles tendons, and gastrocnemius and gastrosoleus (Calf) muscles. Using Dorsi Ramp applies a passive weight bearing stretch to the Achilles tendon and Calf muscles. This type of stretch is often recommended for children suffering from a range of medical conditions where tendon tightness and muscle contracture is affecting the ability to dorsiflex the foot and is possibly leading to an abnormal gait. Examples of conditions where Dorsi Ramp may be applicable include:

- Neurological conditions such as Cerebral Palsy
- Congenital conditions such as Congenital Talipes Equino Varus (Clubfoot)
- Idiopathic conditions, most commonly Idiopathic Toe Walking
- Recovery from trauma and tendonopathy
- Various neuropathies and myopathies

Dorsi Ramp is for use by children able to stand upright and plantigrade (i.e. with their feet flat to the ground) from 18 months of age to not more than 10 years of age and weighing no greater than 63.6kg (140 pounds).

How to Use

To achieve a passive weight bearing stretch to the Achilles tendons and Calf muscles place Dorsi Ramp on a solid level surface, select the correct slope angle and stand upright as directed (1 to 6).

Changing the Slope Angle

Select the most appropriate slope angle, which can be 10, 15 or 20 degrees. To switch between 10 and 20 degree slope angles pick up Dorsi Ramp and turn it upside down. When lifting Dorsi Ramp hold from each end to avoid trapping fingers. The slope angle is clearly indicated. The 15 degree position is achieved by positioning the Dorsi Ramp with the 20 degree slope facing upwards and placing the orange coloured insert on top (1).

Choosing the Correct Slope Angle

Carefully monitor the child to determine the most appropriate angle. Do not use Dorsi Ramp at slope angles that cause discomfort to the child. Excessive stretching can lead to Calf and tendon strains or tears.

A child will be coping well with a slope angle if they can stand on Dorsi Ramp with a normal looking, straight and upright posture for 2 minutes without signs of discomfort; If the child needs to push their bottom out, raise their heels or shows signs of discomfort the slope angle is too high and should be reduced. If the child is not comfortable with the 10 degree slope angle then do not use Dorsi Ramp.

Start with the lowest slope angle and aim to progress to greater angles over a period of weeks, months or even years as appropriate to the user's condition. When the slope angle is increased monitor the child carefully for signs of discomfort; start with short duration stretches before building up to longer ones. There may be days or periods when a child's ability to dorsiflex reduces. This is normal in fast growing children. If the child shows signs of not coping with the slope angle then reduce it.

Frequency and Duration of Stretches

Daily use is recommended. Ideally stretches should last for at least two (2) minutes, but not longer than 30 minutes. Stretches for as little as 30 seconds may still be beneficial. If possible encourage use several times per day.

Correct Foot Position and Posture

The child should stand on the ramp with their heels close to the lowest corner of the slope. The feet can be pointing directly forward or angled outwards as the slope angle is constant in all positions (2).

Should Shoes be Worn When Using Dorsi Ramp?

Children who overpronate or are flatfooted should use Dorsi Ramp with footwear that provides support for the longitudinal foot arch. Ideally such footwear should include custom orthotics. Children who do not overpronate or who tend to supinate (such as children with clubfoot) can use Dorsi Ramp in just socks or bare feet. Shoes, slippers or sandals can always be used if preferred (3).

Encouraging Children to Use Dorsi Ramp

The greatest benefits will be achieved with regular daily use - even for relatively short periods of time.

Persuading young children to engage in regular physiotherapy is not easy. It is best to incorporate the use of Dorsi Ramp into everyday activities and to create incentives and positive associations with its use. Try to avoid situations that lead to negative associations as these will lead to resistance. The best examples of use are when the child is so distracted by what they are doing that they become unaware they are undergoing physiotherapy. Here are some examples:-

Daily Teeth Cleaning: Dorsi Ramp can be placed in front of the bathroom sink and used during the morning and evening teeth cleaning routine. Dorsi Ramp acts as a step, which helps young children comfortably reach the sink. Try using a kitchen timer and insist the child brushes their teeth for 2 minutes each morning and evening (4).

Computer Games: Incentivise children to use Dorsi Ramp by allowing extra time on computer games or other fun activities (5).

Kitchen worktops: Dorsi Ramp acts as a step helping children to comfortably reach the worktop to help out in the kitchen (6).

Storage and Cleaning

Do not store in direct sunlight. Store at temperatures between 5 and 30 degrees Centigrade. Do not stack items on top of the Dorsi Ramp. Plastics will yellow in time, this is normal. Clean the unit using a damp cloth only, do not use bleach, solvents or abrasive cleaning pads. Do not use polish.

Warnings and Precautions

Use only as directed in this guide. Stand only on the sloping surfaces. Do not stand on edges as product may tilt or tip over. If in doubt about how or whether to use Dorsi Ramp please consult a qualified medical practitioner with knowledge of the user's condition. Do not use ramp at slope angles that cause discomfort to the child. Not for use by children less than 18 months of age or unable to stand upright with plantigrade feet. Not to be used by children over 10 years of age or weighing greater than 63.5Kg (140 pounds). Children using Dorsi Ramp must be supervised by an appropriate adult carer at all times. When lifting Dorsi Ramp hold from each end to avoid trapping fingers. Do not drop the Dorsi Ramp. Do not use if damaged or cracked. If in doubt inspect Dorsi Ramp for signs of damage. The Dorsi Ramp is not a toy and should not be used as such. Do not use outdoors.

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Black Robin Designs is a registered trademark of C-Pro Direct Ltd.