

EASY CLICK MD FOOT ABDUCTION BAR

Model EZCLKMD Small and Regular

Indications for Use:

After your child's foot has been corrected, the Easy Click foot abduction brace prevents relapse by holding the foot in the corrected position.

NOTE: Never apply this brace to an uncorrected foot. The brace does not correct clubfoot, it only holds the correction achieved by the Ponseti Method of treatment (a method that uses a series of serial cast to gradually correct clubfoot).

Instructions for Use:

The Easy Click Bar should be worn 23 hours a day for the first 3 months and then at nighttime and naps for 2 to 4 years. Bracing is critical in maintaining the correction of clubfeet. If the brace is not worn as prescribed, there is a near 100 percent recurrence rate.

Warnings and precautions:

Never use Loctite or other superglue-based thread lock on the screws in the Dobbs Bar. The superglue/Loctite will react with the polycarbonate plastic and cause the plastic to weaken and break.

Setting Up the Easy Click Bar

Setting the bar width:

The Easy Click Bar should be set so that the width of the bar is equal to the shoulder width of the child. Measure the shoulder width of the child from the left outside shoulder to the right outside shoulder. Adjust the length of the Easy Click Bar so that the length of the bar is equal to your shoulder width measurement from mid-heel of the left footplate to mid-heel of the right foot plate. It is better to have the bar a little wider than shoulder width rather than too narrow. It is uncomfortable for the child if the bar length is too narrow.

Using the hex key provided in the package, loosen the hex head screws on the center clamp and slide the width of the bars open until they match the shoulder width of the child. Once the width is correct, retighten the hex head screw(s) to lock the bars firmly in place.

Setting the external rotation:

The affected foot or clubfoot should be set at 60 degrees of external rotation, the none affected side should be set at 30 degrees or external rotation. For bilateral clubfoot both feet should be set to 60 degrees. The angle of external rotation should match the degree of external rotation achieved in the last cast. For example, if the physician was only able to correct the affected foot to 50 degrees then the brace should also be set to 50 degrees. *Loosen the hex head screws found on the ends of the bar and rotate the quick clips outward until the pointer on the black Base lines up with the correct degree setting on the bottom of the quick clip. Retighten the hex screws to lock in the desired angle of external rotation.*

Attaching the AFO's

Once the angle of external rotation is set, simply slide the black piece labeled "PUSH" into the recessed channel on the bottom of the AFO. The piece labeled "PUSH" will pop up and lock the AFO in place. To take the AFO off,

*press firmly on the part labeled "PUSH" while pulling backwards on the bar until it slides free of the channel.
Always, put the child's shoes on first then slide them onto the Easy Click Bar.*