

# Product Packaging Medical Device Labelling

Standard format used for most C-Pro Direct products (100mm x 50mm)

**Main label attached to packaging**

GS1  
Data  
Matrix



(01) 5 0562930 10924 6  
 (17) 260203  
 (10) 30222  
 (240) 120120GXXXSLN

GTIN  
 Expiry date (YYMMDD)  
 Lot Number  
 C-Pro product code

Product Description

120STJ 120TTJ GREY XXX-SMALL LEFT  
ADM

Date of manufacture

2022-03-20

Rev: ADMv004

Product Revision



C-Pro Direct Ltd. UK. TN8 6HF

EC REP C-Pro Direct Ireland Ltd. Dublin. D06124

Manufacturer / EU Representative

**Mini label attached to product when larger label attachment is not possible**

120120GXXXSLN C-Pro product code  
 30222 Lot Number  
 260203 Expiry date (YYMMDD)



ADM Ankle Foot Orthosis Label



(01) 5 0562930 05746 2  
(17) 250530  
(10) 210530  
(240) 3GGREVA

GTIN LOT

ADM SANDAL SIZE 3 GREEN AND GREY RIGHT  
EVA

**3 RIGHT**


2021-05-30

C-Pro Direct Ltd. UK. TN8 6HF



C-Pro Direct Ireland Ltd. Dublin. D17PA46

# ADM External Rotation Bar and Bar Clip Labels


Large Bar Kit 1 250mm to 310mm (4 bars)



(01) 5 0562930 14656 2  
(17) 260422  
(10) 290422  
(240) L BK1

GTIN  LOT 

2022-05-03 Rev: v002




CE



C-Pro Direct Ltd. UK. TN8 6HF  
C-Pro Direct Ireland Ltd. Dublin. D17PA46

EC REP


Small Bar Clip Left 60 degree abduction



(01) 5 0562930 06422 4  
(17) 260422  
(10) 290422  
(240) SBCL60

GTIN  LOT 

2022-05-03 Rev: v002



CE

C-Pro Direct Ltd. UK. TN8 6HF  
C-Pro Direct Ireland Ltd. Dublin. D17PA46

EC REP

# D-Bar Product Label



GTR (01) 50562930147934  
LOT (17) 251006  
(10) 100621  
(240) 69MKSA

## DBar Europe 6" to 9" Markell Spring Assist Bar



2021-10-06

Rev: MDv004



C-Pro Direct Ltd. UK. TN86HF



C-Pro Direct Ireland Ltd.  
57 Amiens Street, Dublin 1, D01YY11, Ireland.

# Dorsi Ramp Product Label

This format is used for Dorsi Ramps manufactured prior to the MDR enforcement date

	Manufacturer:	C-Pro Direct Ltd
	Catalogue No:	DORSIRAMP
	Lot number:	CPD004R102
	Manufacture date:	September 2018

C-Pro Direct Ltd, 7a Enterprise Way, Edenbridge,  
Kent. TN8 6HF. UK



5 0 6 0 5 9 1 5 1 2 5 7 8



## Iowa Braces Product Label



C-Pro Direct is a Distributor of the Clubfoot Solutions Iowa Brace. This is an example product label applied to the end of the box

The EU Authorised Representative is Emergo Europe



The LOT number label is applied to the underside of the box

## Despatch Label


This label is attached to Product Packaging and / or Product Cartons at time of despatch to provide customers with order / patient details

**ADM SANDAL SIZE 4 GREEN AND GREY  
LEFT EVA**

*Cust Ref:*

*Order No:* 110326      *Line No:* 1

*Stock Code:* 4GGLEVA



**5056293005753**

C-Pro Direct Ltd, 7a Enterprise Way, Edenbridge Kent,  
UK. TN8 6HF. Enquiries@c-prodirect.com

Product description

Customer's patient ref (if provided)

Sales Order number

Sales Order line number

Product code

GTIN / UDI

Despatch site address (UK address shown)

## Wasp Label (linear)

This label may be attached to Product Packaging and / or Product Cartons to aid internal stock control



**Note:** C-Pro Direct policy is to reduce the number of labels on packaging. When our stock management system software fully supports GS1 Data Matrix Parsing the use of Wasp labels will be reduced.



## Wasp Label (linear or GS1 Data Matrix)

One of these labels may also be attached to Product Packaging and / or Product Cartons to aid internal stock control



The diagram illustrates the structure of a Wasp Label. It consists of three linear barcodes stacked vertically. The top barcode is labeled 'SBK1'. The middle barcode is labeled '130122'. The bottom barcode is labeled '2027-01-13'. Blue arrows point from the labels 'Product code', 'Lot / Batch', and 'Date Code (expiry)' to the respective barcodes. Below the barcodes is a section titled 'Small Bar Kit 170mm to 250mm (5 bars)' which contains a GS1 Data Matrix barcode and a table of data fields.

Length	Value	Field Name
(01)	5056293014649	GTIN
(17)	27.01.13	Date
(10)	130122	LOT
(240)	SBK1	Alt